




U.S. Food and Drug Administration

Enabling Sequence-based Technologies for Micro Diagnostics:



FDA Division of Microbiology Devices Perspective

Peyton Hobson, Ph.D.
Division of Microbiology Devices
OIR/CDRH/FDA/HHS

John.hobson@fda.hhs.gov


Disclaimer

- Sequence-based diagnostic devices for the Microbiology Laboratory are raising new policy / regulatory issues; thoughts presented here are preliminary and do not represent finalized FDA policy
- Pre-submission for outstanding questions

Microbial Dx Sequence Based Technologies

1. Possible approaches to validation studies for High Throughput Sequencing (HTS) systems
 - Metagenomic vs. Targeted (custom amplicon)
 - Use of sequence outputs in combination with database to evaluate performance
 - HTS as a comparator for regulatory submissions
2. Inter-Agency Working Group on Feasibility
3. NIST Microbial Reference Materials

1. Potential Validation Strategies



I. FDA Regulatory Oversight

II. Processes Involved for Comparator Database (Micro DB)


- System approach – collection to result
- Impact of format on validation strategies
 - Metagenomic vs. Targeted (custom amplicon)

System Approach ??

- Sampling
- “pre-analytical” Steps
 - Lysis
 - Nucleic Acid Isolation
 - Library construction
- Amplification
- Sequence Detection/Read
- Assembly
- Database Query/Algorithm
- Database

Evaluation of Diagnostic Devices

FDA's general concept of diagnostic device evaluation



Problem: each possible organism needs confirmation by reference method (ref. positive or negative)

Metagenomic Sequencing Validation Studies

Clinical Evaluation (multisite)

Clinical Performance

- Multisite (end-user)
- Designed to test the Data Analysis Pipeline – Final call determination
- Performance could be established through comparison to database


Precision/Reproducibility

- Panel of representative microorganisms
- Microbial Reference Standards (NIST)

Analytical Studies (in-house)

- Specimen Type and Handling
- Library Prep
- **LOD ***
- **Inclusivity ***
- Interference
 - Specimen
 - Chemistry
- Contamination/Carry-over

HTS as a Comparator



- Cornerstone of all microbiology regulatory submissions
- Accomplished using a comparator device
 - Can be a composite of multiple technologies
 - Burden increases with increase in analytes (especially relevant for multiplexed technologies)
 - May not be uniform

A standardized method that can accurately establish the presence/absence of a microbe in a given sample would be ideal

2. Inter-Agency Working Group on Feasibility

Approach:

- Formed a diverse working group **FDA, NIH-NCBI, NIAID, DTRA, LLNL, and CDC**
- Conducted small pilot study to generate information to evaluate quality of existing sequences in the public domain (In Progress)
 - Identify the pre-existing high-quality deposits, and build from there
 - Will use information to set quality bar for sequence outputs for our ongoing sequencing efforts
 - Utilized existing standards (if available) for technical and isolate metadata –no need to re-invent
 - **Attention given to connecting antimicrobial resistance phenotype to genomic deposits – clinical collection site**

High-level Results Summary


1. **Unbiased benchtop NGS** is currently unable to make confident species calls from mock human samples at clinically-relevant titers
 - Clutter mitigation and/or targeted amplification will be needed**Highly-informative regions were seen at 10 GE/ml using AmpliSeq (LLNL/NMRC)**
2. **Sample carry-over** is a real concern, despite barcoding and aggressive cleaning protocols
Likely to be one of the biggest hurdles for NGS as a clinical diagnostic.
3. **Human sequence** present in some microbial genomes can cause exciting results
LLNL analyzed all human data in NCBI to extract many kmers that are not in the reference human genome -> scalable way of filtering microbial genomes
4. **Screening** against one reference human genome may be insufficient
Some sequencing reagents cause false positive microbial hits if not vigorously masked (phiX174 control, Illumina adapter that has high homology to *E. coli*, NIST controls, etc.)
Aggressive screening of artificial constructs for all reference genome databases.
5. **Restricted reference databases** provide a false sense of confidence about results
Opportunities for false positives, false negatives, and missed detection due to limited taxonomic and/or strain coverage
“Drive-by” NGS analyses are not a long-term solution to the very real problems of genome reference DB curation

Looking ahead: Predictions for Reference Databases

- **Multiple levels of Reference DBs likely**
 - “High quality” genomes only
 - For validation and clinical use
 - “High quality” + other available genomes
 - For testing and development
 - Requires definition of “high quality” that *must* include some draft genomes
- **Extensive screening required**
 - Human and other hosts; chimeras
 - Artificial constructs
- **Separate bacterial, viral, fungal reference DBs**
- **Publicly available (NCBI/EMBL/DDBJ)**

Current Need

Robust, Standardized, and High Quality Microbial Sequence Database in the Public Sector



- Representative Samples
- Metadata
- High quality raw sequences
- Assemblies
- Annotation
- Public Domain

Cover illustration

(Copyright © 2009, American Society for Microbiology. All Rights Reserved.)

3. Microbial Reference Database (MicroDB)

- Identify “gaps” and target sequencing efforts (Funding awarded by FDA/OCET)
 - Raw reads, assemblies, annotations, metadata sent to NCBI and accessible by the PUBLIC
 - Traceable results that could be reevaluated as necessary



Collaborations with Clinical Labs and Repositories

- Children’s National Hospital
- DoD Critical Reagents Program (CRP, USAMRIID)
- FDA-CFSAN, FDA-CBER, FDA-CDER
- DHS National Biodefense Analysis and Countermeasures Center (NBACC)
- The Rockefeller University
- Culture Collections: ATCC, DSMZ

Sequencing Center (UMD IGS)

- Hybrid Approach (PacBio and Illumina)
- Deposit of Raw Reads at NCBI (SRA)
- Deposit of Assemblies at NCBI
- Deposit of Annotations at NCBI
- FDA Interface to Access Data

4. NIST Microbial Reference Materials

- Characterizing Reference Materials for Clinical Applications of Genome Sequencing
- Reference Material (RM) “Material, sufficiently homogeneous and stable with respect to one or more specified properties, which has been established to be fit for its intended use in a measurement process”
- Genomes for NIST RMs
 - **Staphylococcus aureus** – FDA interest
 - **Salmonella typhimurium** – FDA interest
 - **Pseudomonas aeruginosa** – high GC
 - **Clostridium sporogenes** - low GC
- Contact: Nathaniel Olson, Ph.D. and Justin Zook, Ph.D. from NIST

Conclusions

- To realize HTS as a regulated device for microbial diagnostics it is necessary to ensure that comparator information used for interpretation of results is of suitable quality and include appropriate metadata. Efforts are underway...
- There will be difference in the validation strategies for targeted vs. metagenomic sequencing applications
- Will enable a streamlined approach for regulatory evaluation
 - HTS may be another tool for commercial developers to utilize in the evaluation of other types of devices, especially for multiplexed devices
- Highlighted as a need by commercial developers and clinical end-users
 - For clinical adoption really need a simple solution (minimize bioinformatic staffing, etc.)
- FDA Division of Microbiology Devices will continue with current efforts to augment existing sequence information in the public domain



Thank You

FDA Micro Team

Peyton Hobson, Brittany Goldberg, Tamara Feldblyum, Uwe Scherf, Sally Hojvat

Collaborators

IGS@UMD

Lisa Sadzewicz, Luke Tallon, Naomi
Sengamalay

Rockefeller University

Bryan Utter, Douglas Deutsch

Children's National Medical Center

Brittany Goldberg, Joseph Campos

DOD-CRP

Shanmuga Sozhamannan, Mike Smith

DOD-USAMRIID

Tim Minogue

NBACC

Adam Phillippy, Nick Bergman

ATCC

Lisa Kerrigan

DSMZ

Cathrin Sproer

LLNL

Tom Slezak

NIH-NCBI

Bill Klimke, Martin Shumway, David Lipman

NIH-NIAID

Vivien Dugan, Maria Giovani

DTRA

Matt Tobelmann, Chris Detter, Eric

VanGieson, Nels Olsen

CDC

Duncan MacCannell

FDA-CFSAN

Eric Brown, Marc Allard, Andrea Ottessen

NMRC

Ken Frey



Enable regulatory and clinical pathway for sequencing based devices in the Microbiology Laboratory



FDA OIR Micro Team