

REGULATORY FLEXIBILITIES – TESTS FOR RARE BIOMARKERS

MEGAN K. DOYLE, JD, MPH

GLOBAL REGULATORY AND R&D POLICY LEAD, ONCOLOGY, DIAGNOSTICS, DIGITAL HEALTH




DISCLAIMER

The views expressed herein represent those of the author and do not necessarily represent the views or practices of the author's employer or any other party.


This presentation does not constitute legal advice of any type.

CHALLENGES IN DEVELOPING CDX FOR RARE BIOMARKERS



- **Rarity of biomarker or condition often means high unmet need**
- **Rarity makes it difficult to find patients for the trial and samples**
- **FDA expectations have not varied based on rarity of biomarker**

FRIENDS OF CANCER RESEARCH EFFORT




Source: Friends of Cancer Research, Expedited Development of Diagnostics for Therapies Targeting Rare Biomarkers or Indications, https://friendsofcancerresearch.org/wp-content/uploads/Dx_For_Rare_Biomarkers_White_Paper.pdf

WHITE PAPER & KEY THEMES

Propose a framework that would facilitate:

- (1) Enrollment of patients in an efficient manner while maintaining clinical trial integrity
- (2) Approval of a CDx based on requirements that consider the benefit-risk profile and feasibility of obtaining samples



Source: Friends of Cancer Research, Expedited Development of Diagnostics for Therapies Targeting Rare Biomarkers or Indications, https://friendsofcancerresearch.org/wp-content/uploads/Dx_For_Rare_Biomarkers_White_Paper.pdf

PROPOSED MINIMUM REQUIREMENTS MAY IMPROVE PATIENT ACCESS TO TRIALS FOR RARE INDICATIONS OR BIOMARKERS

Performance Characteristic	Minimum Requirement*
Concordance (Sensitivity, specificity, accuracy)	30 biomarker negative samples A range up to 30 biomarker positive [#] samples If possible 6 known positives (confirmed using an orthogonal method)
Limit of Detection	1 known positive* sample in a serial dilution series with at least 3 replicates at each dilution step
Precision	Repeatability across operators, reagent lots, days, instruments using 2 positive samples per variant type, with one at 1.5x LOD and one at 2x LOD
Limit of Blank	5–10 replicates across 2–3 healthy donor samples using the same sample type
*Requirements and number of samples should be guided by the complexity and prevalence of the biomarker being detected	
[#] Can be a contrived sample	

These metrics can also facilitate bridging studies for CDx

Source: Friends of Cancer Research, Expedited Development of Diagnostics for Therapies Targeting Rare Biomarkers or Indications, https://friendsofcancerresearch.org/wp-content/uploads/Dx_For_Rare_Biomarkers_White_Paper.pdf

REGULATORY FLEXIBILITIES – HOW AND WHEN TO APPLY

Recommendations

- Consider a threshold of 10,000 patients likely to have the disease/condition to define rare biomarkers
- Publish benefit-risk assessment in SSED. Rarity may alter benefit/risk.

Table 2. Example Considerations for Benefit-Risk Analysis

Rare Disease/ Variant/ Required Tissue for Validation	Characteristic Qualifying as Rare
NTRK	• Prevalence of Variant (0.32% across solid cancers)
ROS-1	• Prevalence of Variant (1.0% of lung non-small cell lung cancer)
TNBC after progression on primary therapy that metastasizes to the bone	• Tissue type and availability
Fine Needle Aspirate (FNA) in NSCLC	• Paired biopsy and FNAs from the same patient needed for validation

REGULATORY FLEXIBILITIES PROPOSED IN THE PAPER

Analytical Validation

- For approved NGS panels, FDA should leverage pre-existing platform-level validation
- FDA should not require revalidation of variants in the same location or base pair
- FDA should not require validating a variant across all different types of cancers

Clinical Validation

- FDA should rely on the proposed CDx's clinical performance (based on clinical outcome data) rather than on concordance studies to LDTs used for enrollment, which may have unknown/varying performance.
- FDA should consider whether bridging studies could be part of post-market commitments

Other Regulatory Flexibilities

- Allow use of a prespecified modification plan for new indications
- Waive or shift into post-market certain studies (e.g., interfering substances, reproducibility, bridging studies).
- Allow for post-market collection of real-world evidence

FDA should allow sponsors to provide some combination of the following to supplement clinical samples:

- Contrived samples
- Similar tumor types/sample types
- Representative variant validation
- Prior data that demonstrate adequate analytical validity for their assay

Source: Friends of Cancer Research, Expedited Development of Diagnostics for Therapies Targeting Rare Biomarkers or Indications, https://friendsofcancerresearch.org/wp-content/uploads/Dx_For_Rare_Biomarkers_White_Paper.pdf




OTHER SIMILAR ACTIVITIES

Current Challenges in CDx Co-Development:

1. Evolution of Drug Development – rare mutations, biomarkers, speed
2. Use of Multiple CDx & Approved CDx
3. Biomarker-Driven Therapies Outside of Oncology


Proposals Being Developed to Evolve the CDx Paradigm To Address These Challenges

EXAMPLE: REAL-TIME CDX REVIEW PILOT



EX: LIMITED CLINICAL TRIAL ASSAY APPROVAL WITH CONTROLS

Drug approval with PMC - limited testing available with controls in place



SUMMARY

- **Co-development has evolved – rare biomarkers + accelerated programs are more frequent**
- **New CDx technologies available today may allow for regulatory flexibilities to more readily be applied**
- **Stakeholders should continue to work together and explore ideas for addressing challenges**

QUESTIONS?